



## Specification Airshower WDU

Length:	4,30
Width:	3,45
Height:	2,62
Equipment:	
Bodyshowers:	2
Handshowers:	2
Lightning:	3 /230V/ 55W
Blower:	1 /230V/ 1000 W
Evacuation pump:	1 /230V/ 280 W
Soap dispenser:	2
Strecherstand:	1
Showerpallet:	2
Deconstrecher:	1

**Airshower Wet Dry Unit** is a Decontamination Unit designed for Ambulance service, Hospitals and other first responder teams, etc.. The low weight and the units flexibility make it suitable to be transported by helicopter. The unit has two separate department and the complete unit fulfils the international demands and recommendations concerning the philosophy of Decontamination.

As the unit is inflatable it will take only a few minutes to establish the complete unit.

**Technical specification**

Type of heater:	Flowing through.
Width:	570 mm.
Height:	620 mm.
Length:	1670 mm.
Fuel required:	Diesel.
Portable fuel tank:	20 l.
Fuel consumption:	5 kg/h continuous running.
Effect burner:	59 kW
Electric supply:	220 V (GB 110 V ONLY)
Effect consumption during running:	200 W.
Pressure on input water:	2 bar back pressure valve included.
Water thermostat, shower water:	35° C .
Precision temperature shower water:	+/- 2°C.
Max change of temperature:	+ 30°C.

**Normeca** VT  
SUPPLIER OF MOBILE HOSPITALS

P.O. BOX 404, 1471 LOERENSKOG - NORWAY  
Tel.: +47 67 92 76 00 - Fax: +47 67 97 17 66  
E-mail: mobile-hospital@normeca.com

*Agent Norge:*  
**Normeca A/S**

P.O. Box 67, Kjemveien 67, N-1471 Skaarer  
Tel: + 47 67 90 62 50 Fax: + 47 67 97 17 66